



BM Test HbA1c Calibrator

Calibrator for quantitative in vitro determination of HbA1c in blood on BioMajesty systems

Preparation

1. For reconstitution, add exactly 1mL of distilled water and allow the calibrator to stand 10 minutes at room temperature. After 10 minutes, mix the content by gentle inversion. Reconstituted material is stable for 10 days at 2 - 10 °C. Do not freeze. When refrigerated, allow the calibrator to warm to room temperature before use.
2. When leaving the calibrator for 2 hours or longer, mix well before use.
3. Refer to vial label for concentration of Hb ($\mu\text{mol/L}$) and HbA1c ($\mu\text{mol/L}$).

Warnings and precautions

1. The intended use of this calibrator is for **BM test HbA1c** kit by JEOL only. Do not use this product for any other purposes.
2. Use distilled water at around 20 °C for reconstitution. This product may not be reconstituted by cool water. Human blood used in production of this material is tested and proved to be HBsAg, HCV-Ab and HIV-Ab negative. However, it is recommended to wear protective gloves when handling this product and take the same precautions against infections transmitted via human blood as for patient samples.

Package

1 mL x 2 levels x 3 vials

Storage

2 – 10 °C

Stability

The expiration date is stated on the label.

Manufacturer

JEOL Ltd.

2-11-1 Inadaira, Musashimurayama, Tokyo 208-0023, JAPAN