Cat. No. 014-1271-4

CA CAL™ S

Identification of the IVD reagent

CA CAL™ S : Calibrator for LED AUTO Calibration of CA-Systems.

Intended use

The product is applied to calibrate the sensitivity of LEDs used for blood coagulation analyzer.

Principles of the examination method

CA CAL S is a calibrator product used for calibrating the sensitivity of LEDs used for our blood coagulation analyzer.

Components

Polystyrene Latex	0.01%
Sodium azide	0.05%

Warnings and precautions

- For in vitro diagnostic use only. Never use this on the human body.
- If eye or mouth contact occurs by mistake, take emergency measures such as washing the area with plenty of water, see a doctor if needed.
- If get the calibrator on your hand, quickly wash the area with plenty of water.
- The product contains sodium azide(0.05%). Sodium azide can form explosive azides when contacting heavy metals, so dilute the product using plenty of water when disposing of it.
- Never freeze the product.
- Ensure that dust, refuse, or fungi do not get mixed into the product after unsealing it.
- Do not use an outdated product.

Examination procedure

For actual use of the product, refer to the user instruction manual for CA-Systems.
The target value is given on the enclosed lot-dependent table of value.

Please set and use "SAMPLE CUP CONICAL 4 mL (#424-1251-0)" as the reagent container for CA CAL S if the original vial will not be used on the CA-system.

Storage and shelf life of unopened product

Storage and shelf life after first opening

Store the product in a dark place at 2 to 8°C (do not freeze).
The availability expiry date is indicated on both the box and the bottle label.
Once opened, CA CAL S must be used within one day.

Limitations of the examination procedure

The use of this reagent, respectively control or calibrator, is validated on specific analyzers to optimize product performance and meet product specifications. Please refer to Instructions For Use of your analyzer whether the use of this reagent, respectively control or calibrator, is authorised by Sysmex. Sysmex cannot take the responsibility for patient results received from the use of Sysmex reagents, controls or calibrators on unauthorised analyzers. It is the responsibility of the user to validate modifications to these instructions or use of the reagent, control or calibrator on analyzers other than those specified by Sysmex.

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