Fluorocell™ RET

EN

Identification of the IVD reagent

Fluorocell™ RET

Intended use

Fluorocell RET is to be used to label the reticulocytes in diluted blood sample for the assay of reticulocyte count, reticulocyte percent and platelet count in blood with Sysmex automated hematology analyzers.

Principles of the examination method

A sample volume of a whole blood specimen is introduced into the analyzer where a portion of it is automatically diluted into a 1:257 dilution with CELLPACK DFL. Fluorocell RET is then added and the entire dilution is maintained at a constant temperature for a defined time period in order to label the reticulocytes present in the sample. The labeled sample is then introduced into the sheath flow detector where forward scattered light and side fluorescence are measured allowing the reticulocyte count (RET#), the reticulocyte percent (RET%), the RBC count, and the PLT count to be computed. Consult your specific Sysmex analyzer's Instruction for Use for additional information with respect to procedures and principles for reticulocyte analysis.

Components

Polyethylene	0.03 %
Methanol	7.9 %
Ethylene glycol	92.0 %

Warnings and precautions (For in vitro diagnostic use only)

1. Wear gloves and a lab coat for protection. Avoid contact with skin and eyes.
2. In case of skin contact, rinse immediately with plenty of soap and water.
3. In case of contact with eyes, rinse immediately with water or normal saline, occasionally lifting upper and lower lids until no evidence remains. Obtain medical attention.
4. If swallowed, seek medical advice immediately.
5. Do not breathe vapor. In case of accident or you feel unwell, seek medical advice immediately (show the label where possible).

Warning
H226 Flammable liquid and vapour.
H302 Harmful for the environment.
H371 May cause damage to organs.
P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P260 Do not breathe mist/vapours/spray.
P308 + P311 IF exposed or concerned: Call a POISON CENTER/doctor.
P501 Dispose of contents/container in accordance with local and national regulations.

**16CFR1500 (U.S.A.):**

DANGER
POISON
HARMFUL OR FATAL IF SWALLOWED.
VAPOR HARMFUL.
MAY BE FATAL OR CAUSE BLINDNESS IF SWALLOWED.
CANNOT BE MADE NONPOISONOUS.

Examination procedure

1. Put a Fluorocell RET cartridge in the prescribed position and then connect the Fluorocell RET line. Do not refill the Fluorocell RET cartridge with new reagent.
2. Do not remove the IC tag until disposal. All the product information is managed by the IC tag on the label.
3. After setting, reset of the package is not recommended. Removing the reagent cartridge from the analyzer may cause deterioration of the reagent by contamination and opening in a sealing film.
4. Refer to the analyzer's Instruction for Use for further information.

Storage and shelf life of unopened product**Storage and shelf life after first opening**

1. Fluorocell RET has an unopened stability of 12 months from date of manufacture when stored at 2-35 °C in a dark place. See outer package or label on the cartridge for expiration date.
- Please do not use reagent once frozen.

2. Once opened and installed in the instrument, Fluorocell RET is stable for 90 days.

Performance characteristics**Limitations of the examination procedure**

Performance should be within instrument specifications. Refer to the analyzer's instruction for use. Reticulocyte count may be inaccurate due to the influences of abnormal samples. Confirm the reticulocyte count by reference methods if these conditions are indicated. The use of this reagent, respectively control or calibrator, is validated on specific analyzers to optimize product performance and meet product specifications. Please refer to Instructions For Use of your analyzer whether the use of this reagent, respectively control or calibrator, is authorised by Sysmex. Sysmex cannot take the responsibility for patient results received from the use of Sysmex reagents, controls or calibrators on unauthorised analyzers. It is the responsibility of the user to validate modifications to these instructions or use of the reagent, control or calibrator on analyzers other than those specified by Sysmex.

Additional required equipment

Fluorocell RET is intended for only use in conjunction with Sysmex diluent reagent CELLPACK DFL. Consult the analyzer's Instruction for Use for specific reagent requirements.

Primary sample collection, handling and storage

Fluorocell RET is intended for use with blood specimens collected either by venipuncture or micro-sampling by skin puncture. Specimens should be collected in EDTA anticoagulant (EDTA-K2, EDTA-K3 or EDTA-Na2). Micro-sampling specimens can be diluted directly into the diluent without utilization of anticoagulant, or can be collected into micro-collection containers. Note, that the anticoagulant EDTA-Na2 may not dissolve easily in blood, and thus causing fibrin formation or platelet aggregation in some samples. Thorough mixing is required until all dry anticoagulant is dissolved. See in the analyzer's Instruction for Use for further information regarding sample requirements.

Disposal procedures

1. Tightly seal the spout of a cartridge before disposing to prevent residual reagent solution leaks. You may use tape to secure the spout.
2. Disposal procedures should meet requirements of applicable local regulations.

Manufacturer

Sysmex Corporation
1-5-1 Wakinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan

Authorized representative / Distributors

Americas: Sysmex America, Inc.
577 Aptakisic Road, Lincolnshire, IL 60609, U.S.A.
Asia-Pacific: Sysmex Asia Pacific Pte Ltd.
9 Tampines Grande #06-18, Singapore 528735

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Hazardous for the environment (H302).
Very toxic (H371).
Do not drink (P260).
IF exposed or concerned: Call a POISON CENTER/doctor.
Dispose of contents/container in accordance with local and national regulations.

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Caution
Hazardous for the environment (H302).
Very toxic (H371).
Do not drink (P260).
IF exposed or concerned: Call a POISON CENTER/doctor.
Dispose of contents/container in accordance with local and national regulations.

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