REF AW-184-217

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HISCL™ CA19-9 II Calibrator

Intended use

For In Vitro Diagnostic Use
For calibration of HISCL CA19-9 II Assay Kit.

Components

1. HISCL CA19-9 II C0 (CA19-9 0 U/mL) 1 mL x 1
2. HISCL CA19-9 II C1 (CA19-9 10 U/mL) 1 mL x 1
3. HISCL CA19-9 II C2 (CA19-9 40 U/mL) 1 mL x 1
4. HISCL CA19-9 II C3 (CA19-9 100 U/mL) 1 mL x 1
5. HISCL CA19-9 II C4 (CA19-9 250 U/mL) 1 mL x 1
6. HISCL CA19-9 II C5 (CA19-9 1000 U/mL) 1 mL x 1

Warnings and precautions

1. The human CA19-9 used in HISCL CA19-9 II C1～C5 was tested negative for HIV-1 antibody, HIV-2 antibody, HCV antibody, HBs antigen, but it cannot completely rule out risks for infection. Further, this product was not tested for other viruses.
Suppose that this kit is potentially infectious, and handle it carefully.
2. Each Calibrator contains sodium azide. Since sodium azide reacts with lead tubing and copper tubing to generate metal azides which can explode, use a large quantity of water when disposing of it.
3. In case of contact with the eyes, mouth, or hands, perform emergency treatment such as washing with a large quantity of water. If necessary, consult a physician or follow in-hospital infection manual.

Examination procedure

Dispense about 200 µL (5 drops) of each Calibrator.

Perform procedures according to the package insert for the HISCL CA19-9 II Assay Kit, and the analyzer instruction manual.

Storage and shelf life after first opening

Store at 2 - 8°C. The shelf life is 90 days after opening.
Do not freeze.

Manufacturer

KAINOS LABORATORIES, INC.
38-18, Hongo 2-chome, Bunkyo-ku,
Tokyo 113-0033, Japan

Authorized representatives

Asia-Pacific: **Sysmex Asia Pacific Pte Ltd.**
9 Tampines Grande #06-18, Singapore 528735

Traceability of values assigned to calibrators and trueness-control materials

HISCL CA19-9 II C1 - C5 has been adjusted by in-house standard materials.

Date of issue or revision

7/2018

Printed in Japan

REF Catalogue number

IVD In vitro diagnostic medical device

Manufacturer

Consult instructions for use

Temperature limitation

Use by

Batch code

Biological Risks

TSPS4217 (IA830) ²