**[**REF**]** CK-920-648

## SG Calibrator

|                                                                                                                                                                                                                                                          |  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| <div><b>Identification of the IVD reagent</b></div> <div><b>SG Calibrator</b></div>                                                                                                                                                                      |  |
| <div><b>Intended use</b></div> <div>For in vitro diagnostic use only <p>SG Calibrator is intended for use in the calibration of urine specific gravity refractometer of Fully Automated Urine Chemistry Analyzer.</p></div>                              |  |
| <div><b>Principles of the examination method</b></div> <div>SG Calibrator contains sucrose. Since specific gravity is proportional to refractive index, specific gravity of sample is obtained with calibration curve prepared with SG Calibrator.</div> |  |

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |                                           |                           |                                                                        |                           |                                                                        |                           |                                           |                             |                                                                        |                             |                                                                        |                             |  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------|---------------------------|------------------------------------------------------------------------|---------------------------|------------------------------------------------------------------------|---------------------------|-------------------------------------------|-----------------------------|------------------------------------------------------------------------|-----------------------------|------------------------------------------------------------------------|-----------------------------|--|
| <div><b>Components</b></div> <div>Sucrose:</div> <table> <tr> <td>SG Calibrator L (specific gravity, 1.005)</td><td>0.08 w/v<span> </span>%</td></tr> <tr> <td>SG Calibrator M (specific gravity, 1.015 for JSCP, 1.014 for non-JSCP)</td><td>2.41 w/v<span> </span>%</td></tr> <tr> <td>SG Calibrator H (specific gravity, 1.035 for JSCP, 1.031 for non-JSCP)</td><td>7.06 w/v<span> </span>%</td></tr> </table> <div>Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one (3:1)</div> <table> <tr> <td>SG Calibrator L (specific gravity, 1.005)</td><td>&lt; 0.01 w/v<span> </span>%</td></tr> <tr> <td>SG Calibrator M (specific gravity, 1.015 for JSCP, 1.014 for non-JSCP)</td><td>&lt; 0.01 w/v<span> </span>%</td></tr> <tr> <td>SG Calibrator H (specific gravity, 1.035 for JSCP, 1.031 for non-JSCP)</td><td>&lt; 0.01 w/v<span> </span>%</td></tr> </table> <div>Abbreviations: L<span> </span>; Low, M<span> </span>; Medium, H<span> </span>; High</div> | SG Calibrator L (specific gravity, 1.005) | 0.08 w/v <span> </span> % | SG Calibrator M (specific gravity, 1.015 for JSCP, 1.014 for non-JSCP) | 2.41 w/v <span> </span> % | SG Calibrator H (specific gravity, 1.035 for JSCP, 1.031 for non-JSCP) | 7.06 w/v <span> </span> % | SG Calibrator L (specific gravity, 1.005) | < 0.01 w/v <span> </span> % | SG Calibrator M (specific gravity, 1.015 for JSCP, 1.014 for non-JSCP) | < 0.01 w/v <span> </span> % | SG Calibrator H (specific gravity, 1.035 for JSCP, 1.031 for non-JSCP) | < 0.01 w/v <span> </span> % |  |
| SG Calibrator L (specific gravity, 1.005)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | 0.08 w/v <span> </span> %                 |                           |                                                                        |                           |                                                                        |                           |                                           |                             |                                                                        |                             |                                                                        |                             |  |
| SG Calibrator M (specific gravity, 1.015 for JSCP, 1.014 for non-JSCP)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | 2.41 w/v <span> </span> %                 |                           |                                                                        |                           |                                                                        |                           |                                           |                             |                                                                        |                             |                                                                        |                             |  |
| SG Calibrator H (specific gravity, 1.035 for JSCP, 1.031 for non-JSCP)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | 7.06 w/v <span> </span> %                 |                           |                                                                        |                           |                                                                        |                           |                                           |                             |                                                                        |                             |                                                                        |                             |  |
| SG Calibrator L (specific gravity, 1.005)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | < 0.01 w/v <span> </span> %               |                           |                                                                        |                           |                                                                        |                           |                                           |                             |                                                                        |                             |                                                                        |                             |  |
| SG Calibrator M (specific gravity, 1.015 for JSCP, 1.014 for non-JSCP)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | < 0.01 w/v <span> </span> %               |                           |                                                                        |                           |                                                                        |                           |                                           |                             |                                                                        |                             |                                                                        |                             |  |
| SG Calibrator H (specific gravity, 1.035 for JSCP, 1.031 for non-JSCP)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | < 0.01 w/v <span> </span> %               |                           |                                                                        |                           |                                                                        |                           |                                           |                             |                                                                        |                             |                                                                        |                             |  |
| <div>Specific gravity values of SG Calibrator (L, M and H) are assigned according to the nomogram of the Japanese Society of Clinical Pathology<sup>®</sup> (referred to as JSCP) and the nomogram widely used out of Japan (referred to as non-JSCP). As there are various kinds of non-JSCP nomograms, assigned values for the non-JSCP nomogram are reference that manufacturer does not guarantee. There may be a little discrepancy if they are measured using the non-JSCP nomograms.</div>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                                           |                           |                                                                        |                           |                                                                        |                           |                                           |                             |                                                                        |                             |                                                                        |                             |  |

Specific gravity values of SG Calibrator (L, M and H) are assigned according to the nomogram of the Japanese Society of Clinical Pathology<sup>®</sup> (referred to as JSCP) and the nomogram widely used out of Japan (referred to as non-JSCP). As there are various kinds of non-JSCP nomograms, assigned values for the non-JSCP nomogram are reference that manufacturer does not guarantee. There may be a little discrepancy if they are measured using the non-JSCP nomograms.

##### Warnings and precautions

- Before use, please read the analyzer's Instructions for Use and instructions given in the Safety Data Sheet carefully.
- Do not use reagents once frozen.
- Reagents must not be used after its expiry date.
- Be careful not to cut your fingers when removing the flip cap.
- Wear gloves and a lab coat for protection when carrying out work on or with the reagents.
- Avoid direct contact with skin, eyes and mucous membrane and avoid ingestion. In case of contact with skin, eyes or mucous membrane, rinse immediately with plenty of water, and seek medical advice. If swallowed, seek medical advice immediately.
- EUH208

Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1). May produce an allergic reaction.

EUH210

Safety data sheet available on request.

##### Examination procedure

Refer to the Instructions For Use of the dedicated instrument carefully to obtain instrument-specific instructions for the appropriate use of SG Calibrator.

##### Storage and shelf life of unopened product

- SG Calibrator should be stored at 2-10 °C. Keep away from sunlight.
- The expiration date of the unopened product is printed on the box and labels.

##### Storage and shelf life after first opening

Once opened use SG Calibrator within the same day.

##### Interpretation of results

Refer to the Instructions For Use of the dedicated instrument carefully to obtain instrument-specific instructions for the appropriate interpretation of results.

##### Limitations of the examination procedure

Do not use the SG Calibrator as a control with other measurement techniques such as manual refractometer and specific gravity test strip.
The use of this product is validated on specific analyzers to optimize product performance and meet product specifications.
Please refer to the Instructions for Use of your analyzer to confirm that the use of this product is authorized by Sysmex.
Sysmex cannot take the responsibility for patient results received from the use of Sysmex products on unauthorized analyzers.
It is the responsibility of the user to validate modifications to these instructions or use of the product on analyzers other than those specified by Sysmex.

##### Reagent preparation

Remove the product from the refrigerator and allow it to return to room temperature before calibration.

##### Disposal procedures

Disposal procedures should meet requirements of applicable local regulations.

##### Literature references

1) Saito M., et al.: Japanese Journal of Clinical Pathology, 27: 1026-1032, 1979.

|                                   |                                                                                                                                                                         |
|-----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Manufacturer</b>               | <div><div><div><span></span></div><div>Sysmex Corporation</div></div> <div>1-5-1 Wakinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan</div></div>                       |
| <b>Authorized representatives</b> |                                                                                                                                                                         |
| Europe, Middle East and Africa:   | <div><div><div><span></span></div><div><b>[</b>EC<b>]</b> <b>REP</b></div></div> <div>Sysmex Europe GmbH</div> <div>Bornbarch 1, 22848 Norderstedt, Germany</div></div> |
| Asia-Pacific:                     | <div><div><div><span></span></div><div>Sysmex Asia Pacific Pte Ltd.</div></div> <div>9 Tampines Grande #06-18, Singapore 528735</div></div>                             |

|                                                                                                                                                                                                                                                                                                                           |                                 |         |                                 |         |                              |         |  |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|---------|---------------------------------|---------|------------------------------|---------|--|
| <div><b>Product information</b></div> <div>SG Calibrator (USG-100A) including<span> </span>:</div> <table> <tr> <td>SG Calibrator L (yellow liquid)</td><td>5x10 mL</td></tr> <tr> <td>SG Calibrator M (orange liquid)</td><td>5x10 mL</td></tr> <tr> <td>SG Calibrator H (red liquid)</td><td>5x10 mL</td></tr> </table> | SG Calibrator L (yellow liquid) | 5x10 mL | SG Calibrator M (orange liquid) | 5x10 mL | SG Calibrator H (red liquid) | 5x10 mL |  |
| SG Calibrator L (yellow liquid)                                                                                                                                                                                                                                                                                           | 5x10 mL                         |         |                                 |         |                              |         |  |
| SG Calibrator M (orange liquid)                                                                                                                                                                                                                                                                                           | 5x10 mL                         |         |                                 |         |                              |         |  |
| SG Calibrator H (red liquid)                                                                                                                                                                                                                                                                                              | 5x10 mL                         |         |                                 |         |                              |         |  |
| <div><b>Date of issue or revision</b></div> <div>06/2020</div>                                                                                                                                                                                                                                                            | Printed in Japan                |         |                                 |         |                              |         |  |

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