Cat. No.974-0521-6

STROMATOLYSER® -WH

PRODUCT NAME
STROMATOLYSER® -WH

INTENDED PURPOSE (FOR IN VITRO DIAGNOSTIC USE ONLY)

STROMATOLYSER-WH is a reagent that lyses RBC for accurate WBC count determination, WBC tri-modal size distribution analysis and hemoglobin level measurement. The reagent is colorless transparent and contains no cyanide or azide compound. It is intended for use in conjunction with selected Sysmex Automated Hematology Analyzer.

METHODOLOGY AND PRINCIPLE OF THE METHOD

Blood sample collected in EDTA anticoagulant is diluted with **CELLPACK** in a WBC counting container. Then a fixed volume of **STROMATOLYSER-WH** solution (1 volume of **STROMATOLYSER-WH** to 2 volumes of **CELLPACK**) is added automatically to obtain a final dilution of 1:500. The addition of **STROMATOLYSER-WH** lyses the RBC and so the remaining cell stroma is at a level undetectable by the instrument. At the same time, the WBC membrane is preserved and WBC are stabilized at a level detectable by the instrument. They are then counted by the DC method. Hemoglobin is released during RBC lysis, and is converted to the red methemoglobin. A portion of this diluted sample is transferred automatically to the hemoglobin detector where the absorbance of the red pigment is measured to give blood hemoglobin level.

COMPOSITION	
Organic quaternary ammonium salt	8.5 g/L
Sodium chloride	0.6 g/L

WARNINGS AND PRECAUTIONS

Avoid contact with skin and eyes. In case of skin contact, flush the area with water. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. If swallowed, seek medical advice immediately.

STROMATOLYSER-WH is intended for only use with blood samples diluted in the Sysmex diluent **CELLPACK**. Product performance cannot be guaranteed with use of other diluents.

PROCEDURE

- Use the **STROMATOLYSER-WH** at a temperature of 15-30°C. Measuring at a temperature above 30°C or below 15°C may give inaccurate WBC count, WBC tri-model size distribution analysis and hemoglobin level.
- Loosen and remove the cap of the bottle of **STROMATOLYSER-WH** and connect to the instrument.
- Refer to the Operator's Manual of the instrument for detailed information.

STORAGE AND SHELF LIFE AFTER FIRST OPENING

Store **STROMATOLYSER-WH** at 2-35°C. Unopened **STROMATOLYSER-WH** stored at 2-35°C has a product life of 12 months after the date of manufacture. The expiration date is indicated on the container label. Once opened, product stability is 90 days at 2-35°C. **STROMATOLYSER-WH** displaying any signs of contamination or instability, as indicated by cloudiness or color change, should be replaced. Please do not use reagent once frozen.

PERFORMANCE CHARACTERISTICS AND LIMITATIONS OF THE METHOD

When control blood samples (**EIGHTCHECK-3WP**, **EIGHTCHECK L-3WP**, **EIGHTCHECK H-3WP**) are analysed, the WBC count and hemoglobin level should be within the expected ranges. When a normal fresh blood sample is measured by the whole blood mode for 10 consecutive times, the reproducibility (CV %) of the WBC count and hemoglobin level should be below 3.5 % and 1.5 %, respectively. Refer to the Operator's Manual of the instrument for detailed information. WBC counts may be falsely elevated due to the influences of abnormal samples including the following:

1.Nucleated RBC, 2. Cold Agglutinin Disease, 3. Platelet Aggregation, 4. Cryoglobulins. Consult the Operator's Manual for indications of these conditions. Confirm the WBC count by reference eye count methods if these conditions are indicated.

Hemoglobin measurements may be falsely elevated due to the influences of abnormal samples including leukocytosis, lipemia, and abnormal proteins in blood plasma. Confirm the hemoglobin measurement by plasma replacement or plasma blank procedures if these conditions are encountered.

SPECIMEN

STROMATOLYSER-WH is intended for use with blood specimens collected either by venepuncture or micro-sampling by skin puncture. Venepuncture specimens should be collected in EDTA anticoagulant (EDTA-K2, EDTA-K3 or EDTA-Na2). Micro-sampling specimens can be diluted directly into the diluent without utilization of anticoagulant, or can be collected into micro-collection containers with EDTA anticoagulant for dilution at later time.

Note, that the anticoagulant EDTA-Na2 may not dissolve easily in blood, and thus causing fibrin formation or platelet aggregation in some samples. Thorough mixing is required until all dry anticoagulant is dissolved. See in the instrument Operator's Manual for further information regarding sample requirements.

DISPOSAL PROCEDURES

Disposal procedures should meet requirements of applicable local regulations.

MANUFACTURER	
	SYSMEX CORPORATION 1-5-1 Wakinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan

AUTHORISED REPRESENTATIVE / DISTRIBUTORS

Europe:	  SYSMEX EUROPE GMBH Bornbarch 1, 22848 Norderstedt, Germany
North America:	SYSMEX AMERICA, INC. 1 Nelson C. White Parkway, Mundelein, IL 60060, U.S.A.
Japan:	SYSMEX CORPORATION 1-5-1 Wakinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan

PRODUCT INFORMATION	
STROMATOLYSER-WH (SWH-200A)	500 mL x 3 bottles

DATE OF ISSUE OR REVISION	
December 2008	Printed in Japan

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