REF AU-449-468

UC-CONTROL™

Identification of the IVD reagent UC-CONTROL™

Intended use

For in vitro diagnostic use only
UC-CONTROL is used for quality controlling the Sysmex urine chemistry analyzer (UC-3500 and UC-1000).

Principles of the examination method

- UC-CONTROL is a liquid control composed of UC-CONTROL -L and UC-CONTROL -H.
- The UC-CONTROL covers the following analysis parameters. UC-CONTROL contains the ingredients corresponding to the analysis parameters. For details about the measurement principles of the analysis parameters, refer to the analyzer's Instructions for Use and the package insert of the urinalysis test strip.

Analysis parameters: Urobilinogen (URO), Blood (Hemoglobin) (BLD), Bilirubin (BIL), Ketones (KET), Glucose (GLU), Protein (PRO), pH, Nitrite (NIT), Leucocytes (LEU), Creatinine (CRE), Albumin (ALB), Specific gravity (S.G).

Components

Human urine 10-60 %; Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1) <0.01 %; Acetone <5 %

NOTE: The UC-CONTROL is prepared from human urine. Chemical and biochemical substances as well as other constituents derived from human origin are also contained within this control material.

Warnings and precautions

- All human source material used to manufacture UC-CONTROL was non-reactive for antigens to Hepatitis B (HBsAg) and HIV-1 Ag, negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), when tested on single donor basis by the U.S. Food and Drug Administration accepted methods. Because no known test method can assure complete absence of human pathogens, UC-CONTROL should be handled with appropriate precautions.
 - Always wear protective coat, gloves and goggles when carrying out work on or with the UC-CONTROL.
 - Adhere the warnings and precautions for handling described on the product container, package box, package insert, the analyzer's Instructions for Use and use the product correctly. The reliability of the analysis values cannot be guaranteed if the product is used outside the prescribed instructions and test strips.
 - The product must not be used after its expiration date. The product that is not properly stored may fail to maintain adequate performance even though the product has not reached the expiration date.
 - Do not use the UC-CONTROL if bacterial contamination is observed.
 - Do not dilute the UC-CONTROL or do not add foreign materials to the product. This may affect the measurement results.
 - Do not freeze the UC-CONTROL. Do not use a product that is suspected to have been frozen.
 - Do not replenish the UC-CONTROL or mix the residual liquid after use with new product.
 - Do not transfer or store the product in other containers (such as test tubes), otherwise stability is not guaranteed.
 - When a strong shock or force is applied to the container of the product, the container may be damaged.
 - Do not remove the container after use but immediately dispose of it.
 - Do not leave UC-CONTROL with cap removed. This may affect the measurement results.
 - Once test strip is immersed, UC-CONTROL may not be used for quality control of the Fully Automated Urine Chemistry Analyzer (UC-3500) any more.
 - For measurements of the UC-CONTROL, use the dedicated analyzer, reagents, and test strips. For details of the required reagents and test strips, refer to the analyzer's Instructions for Use.
 - After use, please read the Safety Data Sheet carefully.
- Warning
- H317 May cause an allergic skin reaction.
- P261 Avoid breathing spray.
- P280 Wear protective gloves/protective clothing/eye protection/face protection.
- P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
- P362+P364 Take off contaminated clothing and wash it before reuse.
- P501 Dispose of contents/container in accordance with local and national regulations.

Examination procedure

- Remove the UC-CONTROL from the refrigerator and allow it to return to room temperature before use.
- Mix it by gentle inversion to prevent the liquid from bubbling up.
- The following steps depend on each type of analyzer. Carry out measurements according to the applicable instructions. Immediately after measurements, place the cap on to the bottle and tighten, then put the bottle of UC-CONTROL into the package box and return to the refrigerator.

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