[REF] CN-383-000

UF-CALIBRATOR™

Identification of the IVD reagent

UF-CALIBRATOR™

Intended use

UF-CALIBRATOR is a calibration material developed to adjust the sensitivity of the Sysmex Fully Automated Urine Particle Analyzers. UF-CALIBRATOR is only to be used by trained field service personnel. The calibration parameters include:

WBC (/μL), Cond. (mS/cm).

NOTE: Cond. (mS/cm) is not reportable parameter.

Principles of the examination method

UF-CALIBRATOR is used to adjust or verify the sensitivity of Sysmex Urine Particle Analyzers with respect to following conditions:

- When quality control indicates any change in the sensitivity parameters;
- When a major component of the optical units is replaced;
- When a major electric circuit related to the sensitivity is replaced or repaired.

Components

Control particles 0.4 % (W/W)

NOTE: This product contains Latex particles.

Warnings and precautions

- Follow the warnings and precautions for handling described on the reagent container, package box, package insert or the service manual, and use the reagent correctly. The reliability of the analysis values cannot be guaranteed if the reagent is used outside of the written intended use.
- Reagent must not be used after its expiration date.
- Please do not use reagent once frozen.
- Avoid contamination with dust or bacteria after the bottle is opened.
- Remove residual liquid from the cap and the tip of the bottle with the gauze or similar cloth. Fasten the cap tightly immediately after use, otherwise evaporation will increase the density of the reagent.
- If a small amount of the liquid is left in a bottle after repeating several measurements, discard the remainder and open a new bottle. Do not add new calibrator to the opened bottle or mix the remainder with the contents of the new bottle.
- Use UF-CALIBRATOR immediately after mixing. The particles settle at the bottom of the bottle if left to stand for more than 30 seconds, which leads to an uneven particle distribution and thus a measurement error.
- Wear gloves and a lab coat for protection when carrying out work on or with the reagents.
- Avoid direct contact skin, eyes and mucous membrane, and avoid ingestion. In case of skin contact, rinse immediately with plenty of water. In case of contact with eyes, rinse immediately with plenty of water, and seek medical advice. If swallowed, seek medical advice immediately.
- Contamination with black foreign material occurs infrequently, but they don't affect the product performance.

Examination procedure

Please perform measurement in accordance with the calibration procedure, which is described in detail in the service manual of the dedicated instrument.

- Remove a bottle of UF-CALIBRATOR from the refrigerator, and equilibrate to room temperature (10-30 °C) for 20-30 minutes before use. Prepare a container for the reagent.
- Mix a bottle of UF-CALIBRATOR until there is no particle sediment remaining at the bottom, then invert vigorously a further 20 times.
- Immediately (within 10 seconds) after mixing, press the side of the bottle gently to allow 0.6 mL (13 to 18 drops) of the reagent to drop from the tip nozzle into a new container.
- Immediately (within 10 seconds) after distribution, put the container into the prescribed set position and then start the measurement. Refer to the detailed description of the calibrator analysis procedure in the service manual.
- Discard the container after measurement; any container once used in this measurement must not be reused.
- To repeat the measurement, return to step 2.
- When analysis is complete, place the cap on to the bottle and tighten, then put the bottle of UF-CALIBRATOR in to the package box and return to the refrigerator.

Storage and shelf life of unopened product

UF-CALIBRATOR is to be stored at 2-8 °C and should not be exposed to direct sunlight. When handled in this manner, UF-CALIBRATOR is guaranteed stable until the expiration date. The expiration date of the unopened product is printed on the box and on the labels.

Storage and shelf life after first opening

Once opened, product stability is 7 days at 2-8 °C. After opening UF-CALIBRATOR, store in the box to avoid direct contact with the sunlight.

Performance characteristics

Refer to the analyzer's Instructions for Use.

Limitations of the examination procedure

- The intended use of this product with the instrument is limited to those parameters for which assay values are provided.
- Assay values have been established through exclusive use of Sysmex specific reagents, and are valid only with laboratory use of the same Sysmex reagent system.
- Correct results cannot be obtained in any measurement procedure other than calibration procedure of the instrument.

The use of this reagent, respectively control or calibrator, is validated on specific analyzers to optimize product performance and meet product specifications. Please refer to Instructions For Use of your analyzer whether the use of this reagent, respectively control or calibrator, is authorised by Sysmex. Sysmex cannot take the responsibility for patient results received from the use of Sysmex reagents, controls or calibrators on unauthorised analyzers. It is the responsibility of the user to validate modifications to these instructions or use of the reagent, control or calibrator on analyzers other than those specified by Sysmex.

Disposal procedures

UF-CALIBRATOR does not contain any biohazardous substances. Disposal procedures should meet the requirements of applicable local regulations.

Manufacturer

 **Sysmex Corporation**
1-5-1 Wakinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan

Authorized representatives

Europe, Middle East and Africa:

 **Sysmex Europe SE**
Bornbarch 1, 22848 Norderstedt, Germany
Sysmex America, Inc.
577 Aptakisic Road, Lincolnshire, IL 60069, U.S.A.
Asia-Pacific:
 **Sysmex Asia Pacific Pte Ltd.**
9 Tampines Grande #06-18, Singapore 528735

Product information
UF-CALIBRATOR (UCA-100A) 30 mL x 2 bottles
Date of issue or revision
11/2022

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