[REF] BE-740-265

UF-CONTROL™

Identification of the IVD reagent

UF-CONTROL™

Intended use

For in vitro diagnostic use only

UF-CONTROL contains control particles for use in quality control measurement procedure of the Sysmex Fully Automated Urine Particle Analyzers and the Sysmex Fully Automated Urine Particle Digital Imaging Device. The quality control parameters include:

RBC (/ μ L), WBC (/ μ L), EC (/ μ L), CAST (/ μ L), BACT (/ μ L), Cond. (mS/cm).

NOTE: Cond. (mS/cm) is not reportable parameter.

Principles of the examination method

UF-CONTROL contains particles representing red blood cells, white blood cells, epithelial cells, casts and bacteria. When UF-CONTROL is measured by a correctly calibrated instrument in quality control measurement procedure the results should fall within the ranges shown on the assay value table.

Components

	UF-CONTROL -H 0.4 % (W/W)	UF-CONTROL -L 0.1 % (W/W)
Control particles		
NOTE: This product contains Latex particles.		

Warnings and precautions

- Follow the warnings and precautions for handling described on the reagent container, package box, package insert or the analyzer's Instructions for Use, and use the reagent correctly. The reliability of the analysis values cannot be guaranteed if the reagent is used outside of the written intended use.
- Reagent must not be used after its expiration date.
- Please do not use reagent once frozen.
- Avoid contamination with dust or bacteria after the bottle is opened.
- Remove residual liquid from the cap and the tip of the bottle with the gauze or similar cloth. Fasten the cap tightly immediately after use, otherwise evaporation will increase the density of the reagent.
- If a small amount of the liquid is left in a bottle after repeating several measurements, discard the remainder and open a new bottle. Do not add new control to the opened bottle or mix the remainder with the contents of the new bottle.
- Use UF-CONTROL immediately after mixing. The particles settle at the bottom of the bottle if left to stand for more than 30 seconds, which leads to an uneven particle distribution and thus a measurement error.
- Wear gloves and a lab coat for protection when carrying out work on or with the reagents.
- Avoid direct contact skin, eyes and mucous membrane, and avoid ingestion. In case of skin contact, rinse immediately with plenty of water. In case of contact with eyes, rinse immediately with plenty of water, and seek medical advice. If swallowed, seek medical advice immediately.
- Contamination with black foreign material occurs infrequently, but they don't affect the product performance.

Examination procedure

Please perform measurement in accordance with the quality control measurement procedure, which is described in detail in the analyzer's Instructions for Use of the dedicated instrument.

- Remove a bottle of UF-CONTROL from the refrigerator, and equilibrate to room temperature (10-30 °C) for 20-30 minutes before use. Prepare a container for the reagent.
- Mix a bottle of UF-CONTROL until there is no particle sediment remaining at the bottom, then invert vigorously a further 20 times.
- (1) When using the Sysmex Fully Automated Urine Particle Analyzers, immediately (within 10 seconds) after mixing, press the side of the bottle gently to allow 0.6 mL (13 to 18 drops) of the reagent to drop from the tip nozzle into a new container.

(2) When using the Sysmex Fully Automated Urine Particle Digital Imaging Device, immediately (within 10 seconds) after mixing, press the side of the bottle gently to allow 1.0 mL (23 to 28 drops) of the reagent to drop from the tip nozzle into a new container.
- Immediately (within 10 seconds) after distribution, put the container into the prescribed set position and then start the measurement. Refer to the detailed description of the quality control analysis procedure in the analyzer's Instruction for Use.
- Discard the container after measurement; any container once used in this measurement must not be reused.
- To repeat the measurement, return to step 2.
- When analysis is complete, place the cap on to the bottle and tighten, then put the bottle of UF-CONTROL into the package box and return to the refrigerator.

Storage and shelf life of unopened product

UF-CONTROL is to be stored at 2-8 °C and should not be exposed to direct sunlight. When handled in this manner, UF-CONTROL is guaranteed stable until the expiration date. The expiration date of the unopened product is printed on the box and on the labels.

Storage and shelf life after first opening

Once opened, product stability is 30 days at 2-8 °C. After opening UF-CONTROL, store in the box to avoid direct contact with the sunlight.

Calculation of examination results

When UF-CONTROL is measured by a correctly calibrated instrument in quality control measurement procedure, the results shall fall within the ranges shown on the respective assay sheet. Refer to the description on quality control in the analyzer's Instructions for Use of the dedicated instrument for the criteria for evaluating the measurement results.

Performance characteristics

Refer to the analyzer's Instructions for Use.

Limitations of the examination procedure

- The indication values and control limits of UF-CONTROL are those obtained with the specified Sysmex reagents and are valid as long as the same reagents are used.
- Correct results cannot be obtained in any measurement procedure other than quality control measurement procedure of the instrument.

The use of this reagent, respectively control or calibrator, is validated on specific analyzers to optimize product performance and meet product specifications. Please refer to Instructions For Use of your analyzer whether the use of this reagent, respectively control or calibrator, is authorised by Sysmex. Sysmex cannot take the responsibility for patient results received from the use of Sysmex reagents, controls or calibrators on unauthorised analyzers. It is the responsibility of the user to validate modifications to these instructions or use of the reagent, control or calibrator on analyzers other than those specified by Sysmex.

Disposal procedures

UF-CONTROL does not contain any biohazardous substances. Disposal procedures should meet the requirements of applicable local regulations.

Manufacturer

	Sysmex Corporation
	1-5-1 Wakoinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan

Authorized representatives

Europe, Middle East and Africa:

	Sysmex Europe SE
	Bornbarch 1, 22848 Norderstedt, Germany

Americas:	Sysmex America, Inc.
	577 Aptakisic Road, Lincolnshire, IL 60069, U.S.A.

Asia-Pacific:	Sysmex Asia Pacific Pte Ltd.
	9 Tampines Grande #06-18, Singapore 528735

Product information		
UF-CONTROL (UCK-900A) including		
UF-CONTROL -H (UHC-100):	30 mL x 1 bottle	
UF-CONTROL -L (ULC-100):	30 mL x 1 bottle	

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