REF CE-919-553

UF-Fluorocell™ SF

EN

Identification of the IVD reagent

UF-Fluorocell™ SF

Intended use

For in vitro diagnostic use only
UF-Fluorocell SF is used to mark formed elements in urine or in body fluids for determination of red blood cells, casts, etc. counts with the Sysmex fully automated urine particle analyzers.

Principles of the examination method
A sample volume of urine or body fluids is introduced into the analyzer where a portion is automatically diluted 1/4 with the special Sysmex diluent UF-CELLPACK SF and UF-Fluorocell SF. UF-CELLPACK SF reduces the interferences from amorphous salts and mucus and adjusts the condition of the specimen to enable marking of the formed elements. The entire dilution is maintained at a constant temperature for a defined time period in order to allow UF-Fluorocell SF to mark the formed elements in the specimen. The marked specimen is then introduced into the sheath flow detector where light scatter and fluorescent emission are measured, and red blood cells, casts, etc. counts are determined.

Refer to the analyzer's Instructions for Use for the measurement principle.

Components

Polymerine dye 0.05 %
Ethylene glycol 99.9 %

Warnings and precautions

1. Follow the warnings and precautions for handling described on the reagent container, package box, package insert or the analyzer's Instructions for Use, and use the reagent correctly. The reliability of the analysis values cannot be guaranteed if the reagent is used outside of the written intended use.

2. Do not use expired reagents, as the reliability of the analysis values cannot be guaranteed.

3. Do not use a reagent which you suspect has been frozen.

4. When replacing this reagent, do not refill and use the same container.

5. Handle the reagent with care to prevent air bubbles. If air bubbles are present, the performed analysis may be incorrect.

6. Do not remove the RFID (Radio-Frequency Identification Device) tag until disposal. All the product information is managed by the RFID tag on the cartridge.

7. After setting, reset of the package is not recommended. Removing the reagent cartridge from the analyzer may cause deterioration of the reagent by contamination from the already opened sealing film.

8. Always wear protective coat and gloves when carrying out work on or with the reagent.

9. NEVER make contact of this reagent with the human body. Avoid contact with skin and eyes, and avoid ingestion. In case of skin contact, rinse immediately with plenty of water. In case of contact with eyes, rinse immediately with water or normal saline, occasionally lifting upper and lower eyelids until no evidence remains. Obtain medical attention. If swallowed, seek medical advice immediately.

Warning H302 Harmful if swallowed.

P301+P312 IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.

P501 Dispose of contents/container in accordance with local and national regulations.

16CFR1500 (U.S.A.):

WARNING

HARMFUL OR FATAL IF SWALLOWED.

Additional required equipment

UF-Fluorocell SF is intended for use exclusively only with the Sysmex reagents. If other reagents are used, the product performance of Sysmex instruments cannot be guaranteed.

Refer to the analyzer's Instructions for Use for further information.

Examination procedure

Put a UF-Fluorocell SF cartridge in the prescribed position and then connect the UF-Fluorocell SF. When using UF-Fluorocell SF, do not peel off the sealing film.

Refer to the analyzer's Instructions for Use for further information.

Storage and shelf life of unopened product

1. Store unopened UF-Fluorocell SF at 2-35 °C and avoid exposure to sunlight.

2. The expiration date of the unopened product is printed on the box and labels.

Storage and shelf life after first opening

Once opened and installed in the instrument, UF-Fluorocell SF is stable for 90 days.

Performance characteristics

Performance should be within instrument specifications.

Refer to the analyzer's Instructions for Use for performance characteristics.

Limitations of the examination procedure

Always use the reagent within the specified analysis environment for the instrument. Accurate analysis results may not be obtained if analyzed in an environment that does not meet the specified conditions. The results may be inaccurate due to the influences of foreign substances, bacteria growth, etc.

Refer to the analyzer's Instructions for Use for further information.

The use of this reagent, respectively control or calibrator, is validated on specific analyzers to optimize product performance and meet product specifications. Please refer to Instructions For Use of your analyzer whether the use of this reagent, respectively control or calibrator, is authorized by Sysmex. Sysmex cannot take the responsibility for patient results received from the use of Sysmex reagents, controls or calibrators on unauthorized analyzers. It is the responsibility of the user to validate modifications to these instructions or use of the reagent, control or calibrator on analyzers other than those specified by Sysmex.

Primary sample collection, handling and storage

UF-Fluorocell SF is intended for use with human urine and body fluids.

See the reagent's Instructions for Use for further information regarding sample requirements.

Disposal procedures

1. Dispose of the cartridge after removing residual reagent solution and tightly seal the spout of the cartridge before disposing to prevent residual reagent solution leakage. You may use tape to secure the spout.

2. Disposal procedure of cartridge, reagent solution and waste effluents from the instrument should meet the requirements of applicable local regulations.

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