XN CAL™ PF

EN

Identification of the IVD reagent
XN CAL™ PF

Intended use
XN CAL PF is a calibrator for Sysmex X series. It is designed for the calibration and calibration verification of PLT-F.

Principles of the examination method
XN CAL PF is to be used as a calibrator of the Sysmex X series instrument system. Use of stabilized cell preparations for controlling haematology instrumentation is an established procedure. When handled like a patient sample and assayed in the Calibrator Mode of a properly calibrated and functioning instrument, XN CAL PF will provide values indicated on the assay sheet.

Warnings and precautions
Do not inject or ingest. All human source material used to manufacture XN CAL PF was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, XN CAL PF should be handled with appropriate precautions.

Components
XN CAL PF includes stabilized human red blood cells and a platelet component in a preservative medium.


Storage and shelf life of unopened product
XN CAL PF is to be stored in a dark place at 2-8 °C. Avoid freezing this material. When stored properly the unopened product can be used until the expiration date stated on the label on the vial. Storage outside the recommended temperature range causes damage to the product. Do not use damaged material for calibration and calibration verification.

Storage and shelf life after first opening
Opened vials should be used within 4 hours. Afterwards opened vials should be discarded.

Indications of product deterioration
If XN CAL PF fails to perform within expected results as indicated on the assay sheet, there may be a problem with either the calibrator, the reagents or the instrument in use. Proceed as follows:
1. Determine if the instrument system is operating properly and does not require cleaning or maintenance.
2. Check if the reagent system is within the expiration date, if the reagent system is not contaminated, if the reagents are stored properly, and etc.
3. Determine validity of XN CAL PF (i.e. make sure the expiration date, or verify that it has not been frozen).
4. Assay an unopened vial of XN CAL PF (i.e. verify if the opened vial is used over the period of 4 hours).
5. Report any discrepancies to Technical Services of the nearest Sysmex authorized distributor.

Additional required equipment
XN CAL PF is intended only for use on Sysmex X series and its reagents. Acceptable limits of XN CAL PF were established by solely using Sysmex X series and its reagents and are thus comparable only with values likewise established.

Examination procedure
1. Remove the vial from the refrigerator and equilibrate to room temperature (15-30 °C) for 15 minutes before use.
2. Roll each vial between the palms of your hands for 15 seconds.
3. Holding the vial by the ends between the thumb and finger, invert the vial 20 times end-to-end using a very quick turning motion of your wrist during mixing.
4. Analyze the Calibrator reagent in the instrument according to the instructions for use. The pierceable septum in the vial cap allows closed tube analysis.
5. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to instrument analysis.
6. Return to refrigerator (2-8 °C) storage.
Steps 1-6 must be repeated upon removing the sample from the refrigerator for the entire open-vial time period regardless of the method of analysis (open tube, cap piercing, auto sample or manual sample).

Preparation of the calibrator
Since XN CAL PF is containing minimal stabilized cells, utmost care is to be exercised in handling the material.

Calibration procedure
See instrument instructions for use (Chapter 12).

Evaluation of results
An analyzer mean value deviating from the acceptable limits indicates the need for recalibration. However, before changing calibration settings the results of the internal quality control (using Sysmex control bloods) and of the external quality control (SNCS-IQAS ONLINE or external third-party inter-laboratory test) have to be taken into account. If results are contradicting it is to be checked that the calibrator is intact. The justification for recalibration shall be given careful consideration.

Metrological traceability and literature references
The assigned value for XN CAL PF is determined by X series who are calibrated with fresh human blood against the following international conventional reference measurement procedures:

PLT: Platelet counting by the RBC/Platelet Ratio Method – A Reference method, ICSH Expert Panel on Cytometry and ISLH Task Force on Platelet Counting, Am J Clin Pathol 2001; 115: 460-464.
Determined from the RBC/PLT ratio performed by fluorescence flow cytometry with platelets labeled with monoclonal antibodies.

Disposal procedures
This product should not be disposed in general waste but should be disposed with infectious medical waste. Disposal by incineration is recommended.

Manufacturer
 **Sysmex Corporation**
1-5-1 Wakinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan

"Producer" on OEM-Basis: STRECK
7002 S. 109th Street La Vista, NE 68128, U.S.A.

Authorized representative / Distributors
Asia-Pacific: **Sysmex Asia Pacific Pte Ltd.**
9 Tampines Grande #06-18, Singapore 528735

Product information
XN CAL PF 3.0 mL / vial

Date of issue or revision
03/2023 Printed in U.S.A.

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