350521-3

e-CHECK®(XE)

PRODUCT NAME
e-CHECK®(XE)

INTENDED PURPOSE
e-CHECK(XE) is intended to be used as a control for complete blood cell count (CBC), white blood cell differential, reticulocyte and nucleated red blood cell (NRBC) parameters on Sysmex X-SERIES instruments.

METHODOLOGY AND PRINCIPLE OF THE METHOD
e-CHECK(XE) is to be used as a haematology control blood for the quality control of the X-SERIES instrument system. Use of stabilized cell preparations for controlling haematology instrumentation is an established procedure. When handled like a patient sample and assayed in the Quality Control Mode of a properly calibrated and functioning instrument, **e-CHECK(XE)** will provide values within the expected range indicated on the assay sheet.

COMPOSITION
The product **e-CHECK(XE)** includes stabilised human red blood cells, human white blood cells, a platelet and nucleated red blood cell component in a preservative medium.

WARNINGS AND PRECAUTIONS
Do not inject or ingest.
All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.

STORAGE AND SHELF LIFE OF UNOPENED PRODUCT
e-CHECK(XE) is to be stored closed at 2-8°C. When handled in this manner, **e-CHECK(XE)** is guaranteed stable until the expiration date stated on the package and vials.

STORAGE AND SHELF LIFE AFTER FIRST OPENING
Open vials and vials which have been sampled by cap piercing will retain stability for 7 days if stored at 2-8°C after being re-capped.

INDICATIONS OF PRODUCT DETERIORATION
If **e-CHECK(XE)** fails to perform within expected results as indicated on the assay sheet, there may be a problem with either the control blood, the reagents or the instrument in use. Proceed as follows:
1. Determine if the instrument system is operating properly and does not require cleaning or maintenance.
2. Check if the reagent system is within the expiration date, if the reagent system is not contaminated, if the reagent has not been frozen, and etc.
3. Determine validity of **e-CHECK(XE)** (i.e. make sure the expiration date, or verify that it has not been frozen).
4. Assay an unopened vial of **e-CHECK(XE)** to verify if the opened vial is used over the period of 7 days).
5. Report any discrepancies to Technical Services of the nearest Sysmex authorized distributor.

ADDITIONAL SPECIAL EQUIPMENT
e-CHECK(XE) is intended for only use with the following reagents:
diluents: **CELLPACK, STROMATOLYSER-4DL, RET-SEARCH (II) DILUENT**.
sheath reagent: **CELLSHEATH**.
lysing reagents: **STROMATOLYSER-NR, STROMATOLYSER-NR(L), STROMATOLYSER-FB, STROMATOLYSER-IM**.
Hgb lysing reagent: **SULFOLYSER**.
staining reagents: **STROMATOLYSER-4DS, RET-SEARCH (II) DYE, STROMATOLYSER-NR DYE, STROMATOLYSER-NR(S)**.

PROCEDURE
1. Remove a vial of **e-CHECK(XE)** from refrigerator, and equilibrate to room temperature (18-30 °C) for 15 minutes before use.
2. Mix the vial by end-to-end inversion until all red blood cells are completely re-suspended (approximately 20 inversions) before placement in the instrument sample rack.
3. Analyze **e-CHECK(XE)** in the instrument Quality Control Mode according to the X-SERIES Instructions for Use. The pierceable septum in the vial cap allows closed tube analysis.
4. If the vial has been opened for use in the Open Mode of analysis, wipe the threads of both vial and cap with clean lint-free tissue or gauze before replacing cap after sampling.
5. Recap the vial tightly and return to refrigerator (2-8 °C) storage.

PERFORMANCE CHARACTERISTICS AND LIMITATIONS OF THE METHOD
1. Mean assay values for each parameter of **e-CHECK(XE)** are derived from replicate analyses on whole blood calibrated instrumentation. The assay values are obtained using instrument manufacturer's recommended reagents. The values obtained on **e-CHECK(XE)** prior to expiration should be within the expected range. The expected ranges listed represent estimates of inter-laboratory variation for each parameter. Inter-laboratory variation is usually acquired by instrument calibration, maintenance and operating technique. For this reason, the assay values given are guide-numbers useful for control, and are not absolute assays for calibration.
The user must establish values and expected ranges for parameters not listed on the assay sheet. It is recommended that at least 5 consecutive analyses be performed on a properly calibrated instrument to establish the "assay" mean and standard deviation for each level.
The white cell complex has been treated to enhance their stability; therefore, they will not stain to demonstrate typical cell morphology. A microscopic differential analysis of white blood cells cannot be accomplished with **e-CHECK(XE)**.

The intended use of this product with Sysmex X-SERIES instruments is limited to those parameters for which assay values are provided. Values provided by QC analysis by the X-SERIES instruments but not listed on the assay sheet should have their QC target and limit values set to 0 (zero) unless these values are established and accepted by the operating laboratory.
Assay values and limits have been established through exclusive use of Sysmex reagents, and are valid only with laboratory use of the same Sysmex reagent system.
Performance of the control product was established through analysis using the Quality Control Mode of operation of the X-SERIES instruments. Analysis of the product in the clinical laboratory should follow the same process as indicated in the instrument Instructions for Use.
Slight variations in recovered values may occur between Open Mode and Closed Mode of instrument operation (see instrument Instructions for Use). Therefore, if the user chooses to use the product for control of both modes of analysis the results should be filed separately. Combining Open and Closed Mode data may widen the range of SD and CV performance.

DISPOSAL PROCEDURES
This product should not be disposed in general waste but should be disposed with infectious medical waste. Disposal by incineration is recommended.

LITERATURE REFERENCES

1. Henry, J.B. Clinical Diagnostic and Management by Laboratory Methods. Ed.17. W.B. Saunders. Philadelphia, PA, 1984
2. Wintrobe, M.M. "Clinical Hematology", 8th Edition, Lea and Febiger, Philadelphia, 1981.
3. Department of Labor, Occupational Safety and Health Administration. 29 CFR PART 1910. 1030: Occupational Exposure to Bloodborne Pathogens: Final Rule.

MANUFACTURER
 SYSMEX CORPORATION
1-5-1 Wakinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan
"Producer" on OEM-Basis: STRECK, U.S.A.
Omaha, NE 68128, U.S.A.

AUTHORIZED REPRESENTATIVE / DISTRIBUTORS
Europe:   SYSMEX EUROPE GMBH
Bombarch 1, 22848 Norderstedt, Germany
North America: SYSMEX AMERICA, INC.
1 Nelson C. White Parkway, Mundelein, IL 60060, U.S.A.
U.S. Patents 5,432,089;6,200,500;6,221,668;6,399,388;6,403,377;6,406,915
Japan: SYSMEX CORPORATION
1-5-1 Wakinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan

PRODUCT INFORMATION
e-CHECK(XE) 4.5 mL x 8 vials
DATE OF ISSUE OR REVISION
09/2010

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"Produttore" su base OEM: STRECK, U.S.A.
Stampato in U.S.A.