e-CHECK®(XS)

PRODUCT NAME
e-CHECK®(XS)INTENDED USE (FOR IN VITRO DIAGNOSTIC USE ONLY)
e-CHECK(XS) is a haematology control material primarily for intralaboratory quality control of automated, semiautomated and manual procedures that measure components of the blood. Additionally, e-CHECK(XS) can be used in external quality assessment.
Do not use e-CHECK(XS) to calibrate the instrument!SUMMARY and PRINCIPLE
e-CHECK(XS) is to be used as a haematology control blood for the quality control of the XS-Series instrument system. Use of stabilized cell preparations for controlling haematology instrumentation is an established procedure. When handled like a patient sample and assayed in the Quality Control Mode of a properly calibrated and functioning instrument, e-CHECK(XS) will provide values within the expected range indicated on the assay sheet.ACTIVE INGREDIENTS
e-CHECK(XS) contains stabilized human red blood cells, human white blood cells and a platelet component in a preservative medium.

PRECAUTIONS

Do not inject or ingest.
All components in this product from human sources have been tested and found non-reactive for Hepatitis B Surface Antigen, Hepatitis C Virus (HCV) and HIV1/HIV2 antibody using techniques specified by the FDA. However, in keeping with good laboratory practice, this product should be considered potentially capable of transmitting infectious disease.INSTRUCTIONS FOR USE
1. Remove a vial of e-CHECK(XS) from refrigerator, and equilibrate to room temperature (18-30 °C) for 15 minutes before use.
2. Mix the vial by end-to-end inversion until all red blood cells are completely re-suspended (approximately 20 inversions) before placement in the instrument sample rack.
3. Attach the instrument sample rack to the instrument according to the XS-Series Operator's Manual. The pliable cap in the vial cap allows closed tube analysis.
4. If the vial has been opened for use in the Open Mode of analysis, wipe the threads of both vial and cap with clean lint-free tissue or gauze before replacing each sampling.
5. Recap the vial tightly and return to refrigerator (2-8 °C) storage.STORAGE AND STABILITY
e-CHECK(XS) is to be stored upright at 2-8 °C before and after opening. When handled in this manner, e-CHECK(XS) is guaranteed stable until the expiration date stated on the package and vial. Open vials and vials which have been sampled by cap piercing will retain stability for a minimum of 14 days if stored at 2-8 °C.INDICATIONS OF PRODUCT DETERIORATION
If e-CHECK(XS) fails to perform within expected results as indicated on the assay sheet, there may be a problem with either the control blood, the reagents or the instrument in use. Proceed as follows:
1.Determine if the instrument system is operating properly and does not require cleaning or maintenance.
2.Check the system status. If within the expiration date, if the reagent system is not contaminated, if the reagent has not been frozen, and if the vial has not been opened.
3.Determine validity of e-CHECK(XS) (i.e. verify that it has not been frozen).
4.Assay an unopened vial of e-CHECK(XS) (i.e. verify that the opened vial is used over the period of 14 days).
5.Report any discrepancies to Technical Services of the nearest Sysmex authorized distributor.EXCEPTIONS
The mean assay values for each parameter of e-CHECK(XS) are derived from replicate analyses on whole blood calibrated instruments. The assay values are obtained using instruments manufacturer's recommended reagents. The ranges listed on e-CHECK(XS) are to be compared with the expected range. The expected ranges listed represent estimates of inter-laboratory variation for each parameter. Inter-laboratory variation is usually accounted for by instrument calibration, maintenance and operating technique. For this reason, the assay values given are guide-numbers useful for control, and are not absolute assays for calibration. The user must establish values and expected ranges for parameters not listed on the assay sheet. It is recommended that at least 5 consecutive analyses be performed on a properly calibrated instrument to establish the "assay" mean and standard deviation for each level.LIMITATIONS
The white cell components have been treated to enhance their stability; therefore, they will not stain to demonstrate typical cell morphology. A microscopic differential analysis of white blood cells cannot be accomplished with e-CHECK(XS).
The intended use of this product with Sysmex XS-Series Instruments is limited to those parameters for which assay values are provided. Values provided in OC analysis by the XS-Series do not appear on the assay sheet should have their QC target and limit values set to zero (0) unless these values are established and accepted by the operating laboratory.Assay values and limits have been established through exclusive use of Sysmex reagents, and are valid only with laboratory use of the same Sysmex reagent system.
Performance of the control product was established through analysis using the Quality Control Mode of operation of the XS-Series Instruments. Analysis of the product in the clinical laboratory should follow the same process as indicated in the Instrument Operator's Manual.

Slight variations in recovered values may occur between Open Mode and Closed Mode of instrument operation (see Instrument Operator's Manual). Therefore, if the user chooses to use the product for control of both modes of analysis the results should be fed separately. Combining Open and Closed Mode data may widen the range of SD and CV performance.

REAGENTS REQUIRED BUT NOT PROVIDED
e-CHECK(XS) is intended for only use with diluent: CELLPACK,
lysing reagent: STROMATOLYSER-4DL,
staining reagent: STROMATOLYSER-4DS.DISPOSAL PROCEDURES
This product should not be disposed in general waste but should be disposed with infectious medical waste. Disposal by incineration is recommended.MANUFACTURER
 SYSMEX CORPORATION
1-51 Wakino-hama-Kagandori, Chuo-ku, Kobe 651-0073, Japan*Producer* on OEM-Basis: STRECK, U.S.A.
Omaha, NE 68128, U.S.A.AUTHORISED REPRESENTATIVE / DISTRIBUTORS
Europe: SYSMEX EUROPE GMBH
Bornbarch 1, 22848 Norderstedt, Germany
North America: SYSMEX AMERICA, INC.
1 Nelson C. White Parkway, Mundelein, IL 60050, U.S.A.Japan: SYSMEX CORPORATION
1-51 Wakino-hama-Kagandori, Chuo-ku, Kobe 651-0073, JapanPRODUCT INFORMATION
e-CHECK(XS) 1.5 mL x 4 vialsREFERENCES
1.Henry, J.B. Clinical Diagnostic and Management by Laboratory Methods. Ed.17. W.B. Saunders. Philadelphia, PA 1984.
2.Wintrobe, M.M. 'Clinical Hematology', 8th Edition, Lea and Febiger, Philadelphia, 1981.
3.Department of Labor, Occupational Safety and Health Administration. 29 CFR PART 1910. 1030: Occupational Exposure to Bloodborne Pathogens: Final Rule.DATE OF LAST REVISION
January 2006
Printed in U.S.A.

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4.If the vial has been opened for use in the Open Mode of analysis, wipe the threads of both vial and cap with clean lint-free tissue or gauze before replacing each sampling.

5.Recap the vial tightly and return to refrigerator (2-8 °C) storage.

CONSERVATION AND STABILITY
e-CHECK(XS) is to be stored upright at 2-8 °C before and after opening. When handled in this manner, e-CHECK(XS) is guaranteed stable up to the expiration date stated on the package and vial. Open vials and vials which have been sampled by cap piercing will retain stability for a minimum of 14 days if stored at 2-8 °C.INDICATIONS OF PRODUCT DETERIORATION
If e-CHECK(XS) fails to perform within expected results as indicated on the assay sheet, there may be a problem with either the control blood, the reagents or the instrument in use. Proceed as follows:

1.Determine if the instrument system is operating properly and does not require cleaning or maintenance.

2.Check the system status. If within the expiration date, if the reagent system is not contaminated, if the reagent has not been frozen, and if the vial has not been opened.

3.Determine validity of e-CHECK(XS) (i.e. verify that it has not been frozen).

4.Assay an unopened vial of e-CHECK(XS) (i.e. verify that the opened vial is used over the period of 14 days).

5.Report any discrepancies to Technical Services of the nearest Sysmex authorized distributor.

EXCEPTIONS
The mean assay values for each parameter of e-CHECK(XS) are derived from replicate analyses on whole blood calibrated instruments. The assay values are obtained using instruments manufacturer's recommended reagents.

The ranges listed on e-CHECK(XS) are to be compared with the expected range. The expected ranges listed represent estimates of inter-laboratory variation for each parameter. Inter-laboratory variation is usually accounted for by instrument calibration, maintenance and operating technique. For this reason, the assay values given are guide-numbers useful for control, and are not absolute assays for calibration. The user must establish values and expected ranges for parameters not listed on the assay sheet. It is recommended that at least 5 consecutive analyses be performed on a properly calibrated instrument to establish the "assay" mean and standard deviation for each level.

LIMITATIONS
The white cell components have been treated to enhance their stability; therefore, they will not stain to demonstrate typical cell morphology. A microscopic differential analysis of white blood cells cannot be accomplished with e-CHECK(XS).

The intended use of this product with Sysmex XS-Series Instruments is limited to those parameters for which assay values are provided. Values provided in OC analysis by the XS-Series do not appear on the assay sheet should have their QC target and limit values set to zero (0) unless these values are established and accepted by the operating laboratory.

Assay values and limits have been established through exclusive use of Sysmex reagents, and are valid only with laboratory use of the same Sysmex reagent system.
Performance of the control product was established through analysis using the Quality Control Mode of operation of the XS-Series Instruments. Analysis of the product in the clinical laboratory should follow the same process as indicated in the Instrument Operator's Manual.

Slight variations in recovered values may occur between Open Mode and Closed Mode of instrument operation (see Instrument Operator's Manual). Therefore, if the user chooses to use the product for control of both modes of analysis the results should be fed separately. Combining Open and Closed Mode data may widen the range of SD and CV performance.

REAGENTS REQUIRED BUT NOT PROVIDED
e-CHECK(XS) is intended for only use with diluent: CELLPACK,

lysing reagent: STROMATOLYSER-4DL,

staining reagent: STROMATOLYSER-4DS.

DISPOSAL PROCEDURES
This product should not be disposed in general waste but should be disposed with infectious medical waste. Disposal by incineration is recommended.MANUFACTURER
 SYSMEX CORPORATION
1-51 Wakino-hama-Kagandori, Chuo-ku, Kobe 651-0073, Japan*Producer* on OEM-Basis: STRECK, U.S.A.
Omaha, NE 68128, U.S.A.AUTHORISED REPRESENTATIVE / DISTRIBUTORS
Europe: SYSMEX EUROPE GMBH
Bornbarch 1, 22848 Norderstedt, Germany
North America: SYSMEX AMERICA, INC.
1 Nelson C. White Parkway, Mundelein, IL 60050, U.S.A.Japan: SYSMEX CORPORATION
1-51 Wakino-hama-Kagandori, Chuo-ku, Kobe 651-0073, JapanPRODUCT INFORMATION
e-CHECK(XS) 1.5 mL x 4 vialsREFERENCES
1.Henry, J.B. Clinical Diagnostic and Management by Laboratory Methods. Ed.17. W.B. Saunders. Philadelphia, PA 1984.
2.Wintrobe, M.M. 'Clinical Hematology', 8th Edition, Lea and Febiger, Philadelphia, 1981.
3.Department of Labor, Occupational Safety and Health Administration. 29 CFR PART 1910. 1030: Occupational Exposure to Bloodborne Pathogens: Final Rule.DATE OF LAST REVISION
January 2006
Printed in U.S.A.

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Producent på OEM-basis: STRECK, U.S.A.
Trykt i U.S.A.