

Dade® Ci-Trol® 3

I Revision bar indicates update to previous version.

CE0197

Control plasma in the therapeutic range

Intended Use

Dade® Ci-Trol® 3 is used for control of coagulation tests in the upper therapeutic range of oral anticoagulant therapy. With its reduced activities of factors: FII, FVII, FIX and FX, this control is equivalent to plasma from a patient receiving anticoagulant treatment. Dade® Ci-Trol® 3 is preferentially included in each run as an "abnormal plasma" for the control of accuracy and precision. In addition to control values for the prothrombin time (PT), Dade® Ci-Trol® 3 also provides assigned values for the activated partial thromboplastin time (APTT) to enable monitoring reduced activity of the intrinsic coagulation factors as well as the therapeutic range in heparin therapy. The assigned values were established on mechanical and optical coagulometers using Siemens Healthineers reagents. These values were established with reference to an international standard calibrated against pooled citrated plasma.

Reagents

| Reagent | Description | Storage | Stability |
|------------------|--|--|---|
| Dade® Ci-Trol® 3 | Lyophilized reagent containing: • human plasma ^a | 2–8 °C May be used up to the expiry date indicated on the label if stored unopened. | 2–8 °C: reconstituted, 16 hours ^b ; 15–25 °C: reconstituted, 8 hours ^b |

^a produced from pooled, fresh citrated human plasma

^b closed original vial

Warnings and Precautions

For *in-vitro* diagnostic use only.

For laboratory professional use.

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State through your local distribution representative in which the user and/or patient is established.

Safety data sheets (MSDS/SDS) available upon request.



CAUTION! POTENTIAL BIOHAZARD

Each donor or donor unit was tested and found to be negative for human immunodeficiency virus (HIV) 1 and 2, hepatitis B virus (HBV) and hepatitis C virus (HCV) using either tests that are CE marked or FDA approved for this purpose. Because no known test can offer complete assurance of the absence of infectious agents, all human derived products should be handled with appropriate caution.

Dispose of hazardous or biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with all government requirements.

Summary of Safety and Performance (SSP) is available in the European database on medical devices (see Eudamed public website: <https://ec.europa.eu/tools/eudamed>). In case Eudamed is not available, SSP can be delivered by the manufacturer on request.

Preparing Reagents

Add **exactly 1.0 mL of double-distilled water**, then gently tilt the vial and leave to stand for **15 minutes** with the vial closed. Before use, gently mix once more! Dade® Ci-Trol® 3 can be used in the same manner as freshly collected citrated plasma.

Notes on assigned values

Assigned values are calculated with mechanical and optical devices at 37 °C. As a rule, the results obtained (seconds; % of the norm; INR; % activity) should lie within the confidence range indicated in the lot-specific Table of Assigned Values.

However, the results obtain in seconds are especially dependent on the method, the equipment used, the working technique and the specificity of the reagent, and may therefore fall outside the indicated confidence range in exceptional cases.

Procedure

Materials Provided

| REF | Contents |
|--------|---|
| 291072 | Dade® Ci-Trol® 3 Table of lot-specific Assigned Values |
| | 10 × → 1 mL |

Materials Required but not Provided

| Item | Description |
|---|---|
| Coagulation analyzers ^c , such as: | <ul style="list-style-type: none"> Automated Blood Coagulation Analyzer CA-600 series (CA-600 series) AUTOMATED BLOOD COAGULATION ANALYZER CS-2500 (CS-2500 System) AUTOMATED BLOOD COAGULATION ANALYZER CS-5100 (CS-5100 System) Automated Blood Coagulation Analyzer CN-3000/CN-6000 (CN-3000/CN-6000 System) |

^c Availability of analyzers may vary by country.

Please note that the applications on other analyzers can be validated by the instrument manufacturer in accordance with the requirements of the REGULATION (EU) 2017/746 under their responsibility as long as the intended purpose and performance are not modified.

Technical Assistance

For customer support, contact your local technical support provider or distributor.

Definition of Symbols

The following symbols may appear on the product labeling:

| | | | |
|--|---|--|--|
| | Do not reuse | | Use By |
| | Batch Code | | Catalogue Number |
| | Caution | | Manufacturer |
| | Authorized representative in the European Community | | Authorized representative in Switzerland |
| | Contains sufficient for <n> tests | | Biological Risks |
| | <i>In Vitro</i> Diagnostic Medical Device | | Temperature Limitation |
| | Consult instruction for Use | | Non-sterile |
| | CE marking of conformity | | CE marking of conformity with notified body ID number. Notified body ID number can vary. |
| | Contents | | Reconstitution volume |
| | Level | | Keep away from sunlight and heat |
| | Warning | | Danger |
| | Prescription device (US only) | | Device Identification (UDI) barcode |
| | REACH Authorization Number | | |

Legal Information

Ci-Trol and Dade are trademarks of Siemens Healthineers.

Sysmex is a trademark of SYSMEX CORPORATION.

All other trademarks are the property of their respective owners.

© Siemens Healthineers, 2008–2024. All rights reserved.



Siemens Healthcare Diagnostics Products GmbH

Emil-von-Behring-Str. 76

35041 Marburg

Germany

siemens-healthineers.com